

UA-767 Plus

DIGITAL BLOOD PRESSURE MONITOR

Loving Care,
Same BP Model..
Going 10 Years
Strong



- Irregular Heartbeat Indicator
- Clinically Validated
- Blood Classification Indicator
- 1 Touch Measurement
- Smooth and Fast Measurement
- Accurate to plus / minus 3mmHg
- Last reading recall
- Extra large 3-line display for easy reading



FDA
APPROVED

**Clinically
Validated**

CE 0366

IHB

MEM

AC

Clinical Validation of A&D UA767 Plus




Ezperanza I. Cabral, M.D.
Edwin Wenceslao, M.D.
Philippine Society of Hypertension
Pasig City, Metro Manila, Philippines

Conclusion: These data shows that the A&D UA-767Plus and the UA787 devices satisfy the new ESH recommended accuracy levels for both systolic and diastolic pressure. The performance of the devices is reliable and uniform across subgroups of subject with different blood pressure levels.

How to take BP Measurement

1. Place the cuff on the arm (preferably the left arm). Sit quietly during measurement.
2. Press the START. The cuff starts to inflate. It is normal for the cuff to feel very tight. When inflation is complete, deflation starts automatically and the (heart mark) blinks, indicating that the measurement is in progress.
3. When the measurement is complete, the systolic and diastolic pressure readings and pulse rate are displayed. The cuff exhausts the remaining air and deflates completely.



JNC7 Blood Pressure Classification

BP Classification	SBP(mmHg)		DBP(mmHg)
Normal	< 120	and	< 80
Pre-Hypertension	120-139	or	80-89
Stage 1 Hypertension	140-159	or	90-99
Stage 2 Hypertension	> 160	or	> 100