

# seca 206

## Roll-up measuring tape with wall attachment



- Space-saving roll-up mechanism.
- Simple read-out of results through display window.
- Durable tape made of metal.

## seca 206:

When the quality is right,  
so are the measurements.

**M**easuring rods, stadiometers and circumference measuring tapes are relatively simple devices in routine medical work. They still have to fulfill challenging clinical requirements. Whether the devices are for in-patient or mobile use, in large or small spaces, for everyday or infrequent use. They all share the renowned seca quality. A quality you can touch and see – in devices that don't warp even after years of use, in measuring slides that always glide smoothly into position and in the consistently precise measurements that such qualities guarantee.

### The space saving solution with roll-up mechanism.

The roll-up measuring tape seca 206 with wall attachment. Uncomplicated fixing on wall with just one screw. With the display window in the head piece, it's easy to read-off the height. The durable tape is made of metal.

*The wall mounting of  
the measuring tape  
with roll mechanism  
reduces space require-  
ments drastically.*



## seca 206

### Technical Data

- Measuring range:  
0 – 220 cm / 0 – 87"
- Graduation: 1 mm / 1/8"
- Dimensions (WxHxD):  
125 x 125 x 173 mm / 4.9 x 4.9 x 6.8"
- Weight: 202 g / 0.4 lbs

Exclusive Philippine Distributor:

**WellnessPRO**  
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Ave cor Lourdes St.,  
Ortigas Ctr. Pasig City 1605 Philippines  
Mobile no.: +63 922 896 0964  
Tel. no.: +63 2 636 3580, +63 2 637 0930, +63 2 637 0473  
Fax: +63 2 636 7542  
Email: wellnessproinc@gmail.com, jcjuan@wellnessproinc.com  
[www.wellnessproinc.com](http://www.wellnessproinc.com)