

seca 472

Mobile stand for
cable remote displays
of seca scales.



- Stable due to non-slip rubber feet.
- Sturdy, durable construction.

seca 472:

The stable solution for cable remote displays.

A his one is not so easy to bowl over: the stand seca 472 is the transportable solution for holding the display and control element of the digital scales seca 635, seca 657 and seca 675.

Firm hold and easy to transport.

The stand seca 472 is not only a practical alternative to the table support or wall fixture. Equipped with non-slip rubber feet, it also replaces the former model on transport castors and can be positioned everywhere quickly and securely.

Sturdy, durable construction.

The stable stand is made of high grade steel tubing and is thus of durable quality. The display and operating control element is held firmly by a robust clamp. The height and inclination angle of the display element can be adjusted to individual requirements. It can therefore be operated comfortably both from a seated or standing position.



seca 472

Technical Data

- Dimensions (WxHxD):
460 x 1,035 x 460 mm
- Weight: 7.8 kg

Exclusive Philippine Distributor:

WellnessPRO
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Ave cor Lourdes St.,
Ortigas Ctr. Pasig City 1605 Philippines
Mobile no.: +63 922 896 0964
Tel. no.: +63 2 636 3580, +63 2 637 0930, +63 2 637 0473
Fax: +63 2 636 7542
Email: wellnessproinc@gmail.com, jcjuan@wellnessproinc.com
www.wellnessproinc.com