

seca clara 803

The rediscovery
of clarity.



- Two-component rubber coating. For pleasant footing.
- Four stable rubber feet. For a safe standing.
- 28-millimetre LCD figures. For easy reading.
- More than 170 years of experience. For maximum precision.
- From the global market leader for medical weighing and measuring.

Exclusive Philippine Distributor:

WellnessPRO
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Ave cor Lourdes St.,
Ortigas Ctr. Pasig City 1605 Philippines
Mobile no.: +63 922 896 0964
Tel. no.: +63 2 636 3580, +63 2 637 0930, +63 2 637 0473
Fax: +63 2 636 7542
Email: wellnessproinc@gmail.com, jcjuan@wellnessproinc.com
www.wellnessproinc.com

seca[®]
Precision for health

seca clara 803:

The rediscovery of substance.

First take a look at the clear design and then at the finishing: even our simplest scale is complex in its design and robust in its construction. It gives even people who weigh up to 150 kilograms safe footing.



The rediscovery of feeling good.

Go ahead, pick up our starter scale and take a close look. You'll be impressed by its classic design and its surface gleaming with high-quality materials. The high-quality rubber coating consists of a special combination of two components which lends softness and solidity. The material can stand up to all commercially available cleaning and disinfecting products.

Professional can be so private.

Whether black or white, whether for professional or private use, this seca starter scale satisfies the highest demands made by private persons and anyone who appreciates high-quality workmanship. The seca clara 803 has its strengths on the outside and on the inside. Underneath the weighing platform you'll see longitudinal and cross braces set one centimetre apart, giving the scale unparalleled stability.


Lightening the load.

Thanks to its low tare weight of only 1.4 kilograms and the handle built onto the back, the scale is easy to carry, making it suitable for users on the move. Viewing the measurements is just as easy from the unconventionally big display with 28-millimetre-high figures. You see, seca cuts no corners even with its lowest-priced scale.

28-millimetre LCD figures. For easy reading.



Close ribbing for a sure foothold.



seca clara 803

Technical data

- Capacity: 150 kg / 330 lbs / 23 sts
- Division: 100 g / 0.2 lbs
- Dimensions scale (WxHxD): 315 x 35 x 325 mm
- Weight: 1.4 kg
- Functions: Auto-HOLD, automatic switch-off
- Parts to be delivered: 4 AAA batteries, 4 extension feet for carpets