

# seca 876

## Flat scales for mobile use



- Mother/child-function for trouble-free weighing of small children.
- Large leveling base for high stability.
- Easy to transport – scale itself weighs only 4.2 kg.
- Battery power for sites without electrical outlets.

Exclusive Philippine Distributor:

**WellnessPRO**  
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Ave cor Lourdes St.,  
Ortigas Ctr. Pasig City 1605 Philippines  
Mobile no.: +63 922 896 0964  
Tel. no.: +63 2 636 3580, +63 2 637 0930, +63 2 637 0473  
Fax: +63 2 636 7542  
Email: [wellnessproinc@gmail.com](mailto:wellnessproinc@gmail.com), [jcjuan@wellnessproinc.com](mailto:jcjuan@wellnessproinc.com)  
[www.wellnessproinc.com](http://www.wellnessproinc.com)

**seca**<sup>®</sup>  
Precision for health

## seca 876:

### A lightweight for heavyweight use.

In the seca 876 you've found the right solution if you're looking for a lightweight scale (only 4.2 kg) that's stable, compact, easy to clean and equipped with a mother/child-function.



#### Great stability.

The seca 876 is very light at only 4.2 kg, but also very tough. Equipped with a steel frame, the scale also boasts a top-quality ribbed surface, extra-strength two-component ABS plastic and Plexiglas®. Small children and adults up to 250 kg can be weighed accurately. Especially large leveling bases keep the scale stable on any type of flooring.

#### Light but impact resistant.

The combination of steel, ABS plastics and Plexiglas is not only very stable but also especially light, a quality that is important for those interested in mobile use of the seca 876. Just as important for use away from home is the large leveling base which allows the scale to be positioned without any trouble regardless of the surface underneath it. It's easy to step onto the scale, thanks to the flat construction and slip-proof standing surface.

#### Light, but with Mother-Child function.

The scale is easy to operate even with the display and operating elements integrated in the base. Just a light tap of the toe on the scale surface turns on the scale. Infants and small children can be easily weighed while being held by the mother. All she has to do is weigh herself separately first. Then the mother/child-function ascertains the tare of the weight.



*Flat construction and slip-proof surface for a comfortable foothold.*



*Large leveling base on the bottom of the scale for safe placement.*



*Large display for easy reading.*



*Tip-on automatic.*

## seca 876

### Technical Data

- Capacity: 250 kg / 550 lbs
- Graduation: 100g < 150 kg > 200 g / 0.2 lbs < 330 lbs > 0.5 lbs
- Dimensions (WxHxD): 321 x 60 x 356 mm / 12.6 x 2.4 x 14.2"
- Weight of scale: 4.2 kg / 9.3 lbs
- Power: batteries
- Functions: mother/child-function, tip-on automatic switch-on, automatic switch-off
- Optional: stadiometer seca 217 combined with adapter element seca 437, carrying case seca 414 (for combination with stadiometer and adapter element), carrying case seca 421