

seca 954

Digital chair scale with wireless transmission

EMR integrated



- Immediate printout of measurements compared to statistical norms.
- Adjustable footrests fold away, preventing patient from tripping or injuring their ankles.
- BMI function.
- Smooth running wheels.
- At the press of a button measured results can be transmitted wirelessly to a seca 360° wireless digital printer or a PC.

seca 954:

Consideration, comfort and safety.

How can you gently and carefully weigh bedridden patients or those who are unsteady on their feet? The best choice is a seca chair scale, which you can take directly to the patient. That saves you – and the nursing staff – from unnecessary steps. And with the scale's wireless function, transmission of data to a seca 360° wireless digital printer or a PC is quick and easy.

Swiveling arm and foot rests for easy transfer to scale.

Care and comfort were kept in mind when designing the stable model with a circular steel tube frame. The comfortable seat provides a safe position while the swiveling armrests and adjustable footrests make the chair scale more user-friendly. The locking rear wheels ensure that it is safe to sit down and stand up at all times.

Effortless mobility.

The seca 954 has four smooth running pivoting wheels making use possible even in the narrowest of spaces.

Clearly laid out control element.

The positioning of the large LCD display unit behind the chair permits intuitive operation of functions: TARE, HOLD, and BMI. The practical pouch seca 471 attached to the scale frame stores the power adapter.

Precise results for all cases.

As the zero point of the weighing function can change due to external influences – especially in the case of a highly sensitive scale, it is automatically checked and digitally aligned once every second. To ensure that the seca 954 determines the precise weight of restless patients, it is fitted with an effective damping system. An additional function: the BMI value and therefore the nutritional condition of the person weighed can be ascertained automatically at the push of a button. This mobile chair scale, which can be operated by rechargeable battery or power adapter, provides enormous help in rehab, orthopedics, dialysis, nursing homes and pediatrics.



The ergonomically shaped handles let you push or pull the chair scale in any direction.



Paperless and error-free wireless transmission of measurements.

With the seca 360° wireless technology the seca 954 can transmit measurements to the optional seca 360° wireless digital printer. With the network-capable software solutions seca analytics 115 software or seca emr flash 101 and wireless USB adapter seca 456, your PC can receive and analyze measurements and forward them to an Electronic Medical Record (EMR) system. Thus, the seca 954 is EMR integrated and ready to handle digital patient records while meeting all the demands the future will bring. For more information about EMR-integration, visit: www.seca.com

seca 954

Technical Data

- Capacity: 660 lbs / 300 kg
- Graduation: 0.1 lbs < 330 lbs > 0.2 lbs / 50 g < 150 kg > 100 g
- Dimensions (WxHxD): 22.17 x 35.31 x 38.5" / 563 x 897 x 978 mm
- Weight: 54 lbs / 24.5 kg
- Power supply: Power adapter / rechargeable batteries
- Data transmission: seca 360° wireless technology
- Functions: TARE, pre-TARE, HOLD, auto-HOLD, BMI, auto-CLEAR, lbs/kg switch-over, damping, automatic switch-off, SEND/PRINT
- System compatible with: 360° wireless digital printer advanced seca 466, 360° wireless digital printer seca 465, 360° wireless USB adapter seca 456, software seca analytics 115 and seca emr flash 101, stadiometer seca 274 and seca 264
- RS232 compatible

Exclusive Philippine Distributor:

WellnessPRO
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Ave cor Lourdes St.,
Ortigas Ctr. Pasig City 1605 Philippines
Mobile no.: +63 922 896 0964
Tel. no.: +63 2 636 3580, +63 2 637 0930, +63 2 637 0473
Fax: +63 2 636 7542
Email: wellnessproinc@gmail.com, jcjuan@wellnessproinc.com
www.wellnessproinc.com

Toll free +1 800 542 7322 • Fax +1 888 705 7397 • info.us@seca.com • www.seca.com
North America West • seca corp. • 13601 Benson Avenue • Chino, CA 91710 USA
North America East • seca corp. • 7240 Parkway Drive, Suite 120 • Hanover, MD 21076 USA